E - PEN

User Manual  

---

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# Overview

E-Pen™ is an automated, non-surgical microneedling device designed for use by certified aesthetic practitioners and/or trained clinicians. The device incorporates the EL-TE Cartridge™ a sterile microneedling cartridge, and a sleeve for single use only.

## Indications For Use

E-Pen™ is a microneedling device and accessory intended to be used as a treatment to improve the appearance of facial acne scars in adults aged 22 years or older.

## Power & Pen Treatment Settings

1. Ensure pen position is set to “0” before starting a procedure.
2. Power pen on by pressing and holding the On/Off button.
3. Turn the pen dial counterclockwise to increase the cartridge position and clockwise to reduce the cartridge position in 0.1 mm increments.
4. The  symbol indicates current position.
5. Perform a small “test patch” on a patient before a procedure to determine proper treatment. Start at a setting of about 0.2 mm and adjust in approximately 0.1 mm increments until desired effect is achieved.

## How To Use

### 1) Setting

1. Connect a battery to the hand-piece.

### 2) Needle Cartridge Connection & Removal

Step. 1


Step. 2

Step. 3

1. Ensure the pen is powered off and set the pen position to "0" prior to installing the cartridge.
2. Apply single use, non-latex gloves.
3. Open the cartridge by holding it right-side up and pulling back the tyvek cover.
4. Install the EL-TE Cartridge™ by pressing it down and turning clockwise until the cartridge will turn no further.
5. Remove the EL-TE Cartridge™ by turning the tip counterclockwise until it releases.
6. Dispose of used ELTE Cartridge™ via a Sharps container.

\* If the cartridge doesn't correctly fit into the hand-piece, it may malfunction.


\* If an EL-TE Cartridge™ Becomes inadvertently contaminated before or during installation (i.e. dropped on floor, open/broken package, subjected to possible contamination), discard, and obtain new EL-TE Cartridge™.

### 3) Needle Length Adjustment

Adjust needle length by turning the adjustment ring.

## How To Connect Battery

1. Adjust needle length by turning the adjustment ring.

2. Turn the battery clockwise.

## How To Charge Battery

1. Thread the battery securely into the charging unit.
2. Insert the charger adapter into the charging unit and insert the plug into a power outlet.
3. Red light on the charging unit will turn on.
4. Charge 3-6 hours until light turns green.
5. After charging, unthread the battery from the charger and thread into the E-Pen™

# How To Operate

## 1) Power ON/OFF

Connect battery to hand-piece and turn on the power button for operation.

Press the power button to turn on/off the power.

## Procedure Instructions

1. Have patient review and complete consent form.
2. Explain and answer any questions about E-Pen™ procedure.
3. Clean patient's skin with soap and water, alcohol pads or a cleanser to remove makeup, sunscreen and surface oils.
4. Take "before" pictures of the procedure area.
5. Apply numbing cream to the procedure area. Leave on per the product's instructed time.
6. Cleanse again, thoroughly removing all numbing cream. Wipe skin with alcohol to ensure all numbing agent is removed.

**IMPORTANT:** Anything remaining on the skin may increase the risk of complications post-procedure.

7. Apply the disposable EL-TE Cartridge™ to pen. Sleeve is non-sterile and single use only.

**1**  
Hold Sleeve and insert E-Pen™ body (without cartridge) between the white tab and paper backing.

**2**  
Push E-Pen™ through open end to expose the first ridge (2 mm) of the tip.

**3**  
Peel back the protective top cover and the paper backing.

**4**  
Pull sleeve snug, remove adhesive backing and seal end.

**5**  
Insert new microneedle cartridge into E-Pen™. E-Pen™ is now protected and ready for use.

8. Install the EL-TE Cartridge™. See Step 7 for installation instructions
9. **IMPORTANT:** The EL-TE Cartridge™ is not intended to be filled with liquids or any substance.
10. Continually apply sterile saline as a lubricating topical to the treatment area as needed during treatment session.
11. Select pen position according to approximate thickness of epidermis. See "Power & Pen Treatment Settings". Microneedling should not be used on areas with thin skin, such as the eyelids and periorbital area.
12. For facial procedures, divide into four quadrants. Start with the right cheek, move to the chin/lip/nose, then to left cheek, and finish with forehead.
13. Hold patient's skin taut in each area. Repeat the pattern if the erythema (redness) endpoint is not reached.
14. Observe the effect on the skin and adjust cartridge position accordingly.
15. An optional third pass may be made using gentle, one-directional circular motions in small targeted areas.

## Post-Procedure Instructions


1. Gently use sterile gauze to pat down the affected area.
2. Apply a hydrating balm.
3. It is recommended to avoid other facial aesthetic treatments following the E-Pen™ treatment.
4. Advise patient to avoid sweaty exercise and sun exposure for 72 hours post-procedure.
5. Schedule next appointment after at least 4 weeks.
6. Take "after" pictures before next appointment.

Hold the E-Pen™ perpendicular to the floor, or with the cartridge attachment tip pointing downwards. Use one hand to remove the cartridge and dispose of the cartridge in a sharps container.

Continue to hold the E-Pen™ perpendicular to the floor with the cartridge tip pointed downwards, and pull apart the adhesive strip of the sleeve.

Remove the sleeve by carefully sliding it down the E-Pen™ to prevent soiling the handpiece.

Dispose of the sleeve in a biohazard container. Sleeves are not intended to be reused.

\* After removal of the sleeve users' gloves should be removed, hands cleaned, and a new pair of clean gloves worn before proceeding to the next patient.

**IMPORTANT:** Unexpected complications may occur when products not proven safe for use with microneedling are applied post-procedure.

**NOTE:** Soiled gloves should always be disposed of in a biohazard container. Do not reuse disposable gloves.

**NOTE:** The purpose of a sheath is to provide a covering that helps prevent the transmission of pathogens from one patient to another. E-Pen™ is intended to be used only with provided sleeve.

## Contraindications

Treatment with this device is contra-indicated for patients with any of the following conditions:

- Keloid scars
- History of eczema, psoriasis and other chronic conditions
- History of actinic (solar) keratosis, diabetes, raised moles or warts on targeted area
- History of hemophilia, irregular blood pressure, tuberculosis, liver function issues
- Susceptibility to capillary ectasia due to steroid use for extended periods
- Scleroderma
- Collagen vascular disease
- Cardiac abnormalities, pacemaker, blood clotting problems
- Blood thinning medication
- Active bacterial or fungal infection
- Immunosuppression
- Facial melanosis
- Malignant tumors
- History of any type of cancer or suspicious lesions or moles in treatment area
- Pregnant or nursing women
- Any other medical condition contraindicated by the treating physician
- Have open wounds, sores, or irritated skin in the treatment area(s)
- Have an allergy to stainless steel or anesthetics
- Have a hemorrhagic (bleeding) disorder or hemostatic (bleeding) dysfunction

**NOTE:** This product is not intended for transdermal (under the skin) delivery of topical products such as cosmetics, drugs, or biologics.

## Precaution Before Use

1. Discontinue auto-immune therapies or retinoid (Retin-A) and/or any form of skin treatment 24 hours prior to procedure, under care and direction of a physician.
2. Do not treat active acne, rosacea, or other inflammatory skin conditions.
3. Do not treat subject with piercings in treatment area or open wounds.
4. Not to be administered for 6 months after isotretinoin (Accutane™) regime.
5. Patients with facial outbreaks such as herpes simplex virus, medication must be taken per doctor's instructions.
6. Patients with metal allergies or skin allergies.
7. Subjects on any medications that would affect the characteristics of the skin should be stopped for two months prior to treatment under the care and direction of a physician.
8. Subjects who had a face lift or eyelid surgery within the past year, or dermabrasion, remodeling, deep chemical peels, or any surgical procedure on the treatment area within the previous 3 months.
9. Subjects who had Botox™, collagen, fat, or other methods of augmentation with injected materials in the treatment area in previous 6 months.
10. Universal precautions are necessary during microneedling. Microneedling should not be used within the orbital rim of the eye, such as the eyelids.
11. Subjects who have excessively tanned or sunburned skin from the sun, tanning beds, or tanning creams within previous 2 weeks.
12. Patients are to avoid sweaty exercise and sun exposure for 72 hours post-procedure.
13. Fitzpatrick skin types V-VI, pigment may darken prior to lightening.
14. Patients could experience redness between one to three days.
15. Patients may experience inflammation, itching and burning.

## Cleaning & Storage Instructions

### Following the Procedure:

1. Immerse the device in Cidex® OPA at room temperature for 90 minutes. (Only immerse from the tip of the device to the top of cone).
2. Rinse the cone end of the device (minimum of 3 times) with deionized water.

1

TYPE B APPLIED PART

2

Class 2 Equipment

- |    |   |  |    |   |  |
|----|---|--|----|---|--|
| 3  |   | Alternating current                            | 4  |   | Stand-by                                     |
| 5  |  | Do not throw away with general household waste | 6  |  | General warning sign                         |
| 7  |  | General mandatory action sign                  | 8  |  | General prohibition sign                     |
| 9  |  | Keep away from rain                            | 10 |  | Fragile; handle with care                    |
| 11 |  | This way up                                    | 12 |  | Handle with care                             |
| 13 |  | Manufacturer                                   | 14 |  | Serial number                                |
| 15 |  | Operating instruction                          | 16 |  | Temperature limits 0 °C min. to + 36 °C max. |
| 17 |  | Humidity limit 30% to 85% R.H.                 | 18 |  | Air pressure limit 800 hPa to 1060 hPa       |

## Specifications

Product name	E-PEN™
Model number	DOA-EPEN
Manufactured for	Salient Medical Solutions Inc.
Weight after packaging	1.5 lbs, 1 set/box
Electrical Requirements	AC Adapter In: 100-240Vac, 50-60Hz 0.3A AC Adapter Out: 5.0Vdc; 2.0A
Charge time	3-6 Hours
Working time	4 hours per battery
Speed	7,000 RPM (puncture rate)
Needles	12 pin, 32 g, Stainless Steel, 2.0 mm Maximum Penetration Depth
Operation	Cordless

# Product Warranty

Salient warrants to the original purchaser that the merchandise is free from defects in material and workmanship under normal use and service for a period of one (1) year from date of purchase. The obligation of Salient under this warranty is limited, in its exclusive option, to repair or replace parts and materials which prove to be defective.

One or more of the following events will serve to void the warranty:

1. Defects due to negligence, alteration, or modification, by anyone other than factory authorized personnel.
2. Abuse or misuse.
3. Attempted or actual dismantling, disassembling, service or repair not specifically authorized by Salient.

All merchandise should be inspected for obvious damage upon arrival. If merchandise has been damaged in transit, Salient must be notified within 48 hours. All claims for nonconforming or defective product must be made in writing within 10 days after delivery to the purchaser, and any claims not made within that period shall be deemed waived and released.

IN NO EVENT SHALL SALIENT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES DUE TO ANY CAUSE WHATSOEVER.

No suit or action shall be brought against Salient more than one year after the related cause of action has occurred.

THE FOREGOING CONSTITUTES SALIENT'S SOLE LIABILITY AND THE PURCHASER'S SOLE REMEDY WITH RESPECT TO PRODUCTS SOLD BY SALIENT EXCEPT AS THUS PROVIDED, SALIENT DISCLAIMS ALL WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

## WARNINGS

Electromagnetic waves discharged from mobile phone, radiotelegraphs, wireless toys can cause device malfunction. Keep away any device that can influence this product. Do not turn the device OFF while the device is in direct contact with the patient. Only use saline directly before and after using the device. Never use unauthorized parts. All parts should be replaced by manufacturer's authorized service personnel. Never use flammable cleaning agents (e.g. benzene, paint thinner, etc.).

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All product names, logos, brands, trademarks and registered trademarks are property of their respective owners. All company, product and service names used in this manual are for identification purposes only. Use of these names, trademarks and brands does not imply endorsement.